

hemochroma^{PLUS}

The hemochroma PLUS[®] Analyzer provides an accurate hemoglobin and hematocrit measurement.



- Reagent-free microcuvettes
- Microcuvette shelf life of up to 2 years, even after the canister has been opened
- Cleared for use on patients 6 months or older
- Measurement range 5.0 to 25.6 g/dL
- Small 15 µl sample volume
- 3.25" Color LCD display
- CLIA Waived





Reimbursement	
Test	Hemoglobin
CPT Code	85018-QW
Blood count; hemoglobin (Hgb)	

Simple 3 Step Procedure



• Performance



• Parameters

Specifications	hemochroma+
Sample Type	Capillary Blood/Venous Blood
Reaction Time	3 seconds
Measurement Range	5.0 to 25.6 g/dL
Sample Volume	15 µl
Microcuvettes Stability	24 Months
Storage Condition	15 - 35°C / 59 - 95°F
Tests/box	200

• Analyzer

Specifications	hemochroma+
Dimensions	4.19 x 5.94 x 1.51 in
Weight	0.5 lbs.
Interface	3.25-inch screen External thermal printer
Memory	1,000 results and QC results
Power	4 AA Batteries*

• Order Information

Products	Catalog Number
Analyzer	100-101
Microcuvettes	100-202
Controls (Level 1,2,3)	100-303
Optical System Check	100-404
Analyzer Promo**	100-PROMO

* Batteries and USB Power Cable Included.

**Promo Kit Includes 1 Analyzer, 1 Optical System Check, 100 Microcuvettes.

PS-1461 Rev04 02/2019



Visit us on the web
at:
www.IMMUNOSTICS.com

38 Industrial Way East Ste.1, Eatontown, NJ 07724 USA
sales@immunostics.com
Tel 732-918-0770 Toll Free 800-722-7505 Fax 732-918-0618